

# iPhone 12 and MagSafe Charging Accessories

The iPhone 12 and Magnetic Charging devices can seriously compromise pacemakers and defibrillators. Please physically keep a safe distance between your iPhone 12 and your charging devices while in either Mercy Hospital location.

- Maintain a distance of at least 15 cm/6 inches between iPhone 12 models and MagSafe accessories and medical devices such as implanted pacemakers and defibrillators.
- When wirelessly charging iPhone 12 models with MagSafe accessories, the distance needs to be at least 30 cm/12 inches from medical devices.



If you have any questions, please connect with a member of our care team. Thank you!

For more information, please visit:

<https://support.apple.com/en-us/HT211900>

